Pipettierhelfer / Pipette Controller

maneus[®]

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1. Safety Instructions

This instrument may sometimes be used with hazardous materials, operations and equipment. This Manual does not purport to address every safety issue which may arise during use. It is the responsibility of the user of this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Every user must read and understand this operating manual before operation.
2. Observe general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves.
3. Carefully observe the operating manual and information from reagent manufacturers.
4. Use the instrument only for pipetting liquids, with strict regard to the defined limitations of use. Observe operating exclusions (s. page 21)! If in doubt, contact the manufacturer or supplier.
5. Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Only use suitable vessels.
6. Never use force on the instrument!
7. Use only original manufacturer's accessories and spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
8. Always check the instrument for visible damage before use. If there is a sign of a potential malfunction, immediately stop pipetting. Consult the 'Troubleshooting' section of this manual (see page 32), and contact the manufacturer if needed.

2. Application

The instrument serves to facilitate the pipetting of liquids with graduated pipettes, volumetric pipettes and blow-out pipettes made of glass or plastic within the volume range of 0.1 ml to 200 ml with a suction tube outer diameter < 9.2 mm for measuring liquids. If the instrument is correctly used, the liquid will not come in contact with the instrument.

3. Limitations of use

The instrument is designed for pipetting liquids, observing the following physical limits:

- +10 °C to +40 °C (50 °F to 104 °F)
(of instrument and reagent)
- Vapor pressure up to max. 500 mbar.
Aspirate slowly above 300 mbar, in order to prevent the liquid from boiling.
- Density up to 2.4 g/cm³

4. Operating Exclusions


The instrument must never be used with liquids that may release vapors which attack PP, PTFE or silicone. The instrument is not designed for Pasteur pipettes.

5. Storage conditions

Store the instrument and accessories only in cleaned condition in a cool and dry place.

Storage temperature: -20 °C to +50 °C (-4 °F to +122 °F).

6. Components


1. suction bellows
2. blow-out button
3. valve unit
- 3a. filteradapter
4. pipetting knob
5. membrane filter
6. adapter
7. adapter support

7. Pipetting

7.1 Compress the suction bellows

- Before attaching the pipette, squeeze the suction bellows (1).

The instrument is now ready to aspirate liquid.


7.2 Attaching the pipet

Attention!

Never use force. Glass and plastic pipets can break and cause injury.

- Hold the pipet near the top and insert into the adapter (6) until it fits tightly.
- Once the pipette has been attached, always hold the instrument in a vertical position, tip pointing downwards.


7.3 Filling the pipet

- Immerse the tip of the pipet into the liquid.
- Move the pipetting lever (4) upwards.
- Fill the pipet to just above the desired mark.

Attention!

Do not over aspirate and allow any liquid to enter the instrument.

This would impair the filtering function (5) and reduce the suction capacity. If liquid does enter the instrument, exchange the membrane filter (5) (pages 27-29).

Note:

- The fill/delivery rate is continuously variable with the upward/downward pressure on the pipetting lever.
- In the case of large pipettes > 50 ml the vacuum contained in the suction bellows (1) is not sufficient to draw in all the liquid at once. Therefore, squeeze the suction bellows (1) again and continue drawing up liquid.

7.4 Adjusting the meniscus


- Use a suitable tissue to wipe the pipet tip. Push down the pipetting lever (4) slowly until the meniscus is exactly level with the desired mark.

7.5 Discharging the pipet

Hold the collecting vessel at an angle. Place the pipet tip against the inner wall of the vessel, push down the pipetting lever (4) to dispense the liquid.


- **Free discharge for pipets with waiting time (Imprint e.g. "Ex+ 15 s")**

- As soon as the meniscus stops moving, start the necessary waiting time (e.g. 15 sec).
- Wipe the pipet tip along the wall of the vessel by sliding upwards a few millimeters.

- **Blow-out pipettes (Imprint "ausblasen – blow out"):**

- As soon as the meniscus in the pipette tip comes to a standstill, press the blow-out bellows (2) once.
- Wipe the pipette tip a few millimeters upward along the wall of the vessel.


7.6 After pipetting

To remove the pipet, gently grasp near the top and carefully remove with a twisting motion.

8. Maintenance & Cleaning

With normal use, the instrument is maintenance-free. The exterior may be wiped with a damp cloth. The complete instrument is autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi) with a holding time of at least 15 minutes according to DIN EN 285.

For autoclaving remove suction (1) below.

9. Replacing the membrane filter/ Cleaning the adapter

Attention!


When cleaning and/or exchanging the membrane filter (5), wear hand and eye protection to avoid the risk of splashes and contamination. Remaining liquid can run out.

If liquid has entered the adapter (6) or rather the membranefilter (5) or if the aspirating efficiency has decreased, proceed as follows:

1. Hold the instrument over a suitable vessel to allow all liquid to drain from the pipette.
2. Grip the pipet as near the top as possible and carefully pull it out while turning.
3. Unscrew the adapter support (7).
4. Pull out the membrane filter (5).
5. Remove the adapter from the adapter support (7) by pushing it upwards.
6. To clean adapter (6), carefully rinse it by using an appropriate compatible solvent (e.g., water) in a wash bottle. Blow out the liquid and allow to dry completely.
7. Place the new membrane filter (5) with the thicker end facing downwards into the adapter (6) (with the smaller conical side up).


8. Screw in the adapter support (7) hand tight.
9. The entire instrument incl. membrane filter is autoclavable at 121 °C (2 bar) according to DIN EN 285.

For autoclaving remove suction (1) bellow.

Note:

Not every membrane filter is suitable. Use only genuine manufacturer's recommended accessories. Membrane filters of a pore size of 3 µm are a standard equipment of the instrument. Up to 5 times, the membrane filter can withstand steam sterilization at 121 °C (2 bar absolute, i.e. 1 bar above atmospheric pressure) for 20 minutes. It is the user's responsibility to ensure effective autoclaving.

**Attention!**

Never use force! While mounting the housing it has to be paid attention for that the pipetting lever is loosely hanging. Then place the pipetting lever scheduledly between the valves and close the housing (s. figure). Does the pipetting lever not hang loosely, but has twisted itself, the valve unit can be damaged sustainably!

10. Check for leaks.


10. Replace the valve unit

Attention!

When cleaning and/or exchanging the valve unit, wear hand and eye protection to avoid the risk of splashes and contamination. Remaining liquid can run out.


If liquid has entered the valve unit proceed as described below:

1. Unscrew the suction (1) bellows and the adapter support (7).
2. Take off the upper part of the housing. Carefully pull out the valve system and replace it.
3. Reassemble the instrument in reverse order and carry out leak test.

11. Checking for leaks

When the pipet has been filled, no liquid should drip unless the pipetting lever is depressed. If leaking occurs, carefully pull out the pipet while turning and dismantle the instrument again. If necessary replace membrane filter (5), clean adapter (6) or replace valve unit (3).

12. Troubleshooting

Trouble	Possible Cause	Action to be taken
Suction capacity impaired	<ul style="list-style-type: none">- membrane filter (5) dirty or- valve unit (3) damaged	<ul style="list-style-type: none">- replace membrane filter (5)- replace valve unit (3)
Pipette drips	<ul style="list-style-type: none">- pipette not inserted far enough or- membrane filter (5) not properly positioned or- valve unit (3) damaged	<ul style="list-style-type: none">- carefully push pipet further into securing sleeve or- insert membrane filter (5) or- replace valve unit (3)
Pipette not held properly	<ul style="list-style-type: none">- adapter (6) damaged	<ul style="list-style-type: none">- replace adapter (6)
Compressed suction bellows (1) take in air automatically	<ul style="list-style-type: none">- suction (1) bellows or- valve unit (3) damaged	<ul style="list-style-type: none">- replace suction bellows (1) or- replace the valve unit (3)

Note: If recommendations in the Troubleshooting do not resolve current problems, contact your supplier.

13. Product Markings

Markings	
	General warning sign
	Observe the operating manual
XXZ	Batch number

14. Ordering Data

maneus® 1630500
complete, with membrane filter

Spare Parts:

Spare membrane filter (5), 1670650
3 µm, non-sterile, hydrophobic, PTFE (10 per pack)

Suction bellows (1), 1670201
complet, red

Adapter (6), 1670301
silicone

Valve unit (3) 1655590
complete

Adapter housing (7) 1661281
PP

15. Repair service

Attention!

Only cleaned and decontaminated instruments can be evaluated/repaired.

Clean the instrument thoroughly. If the instrument has been contaminated with infectious, genetically modified or biologically active or radioactive substances, the instrument must be decontaminated in addition.

Return instrument along with a description of the malfunction and dispensed liquids. Return instrument along with a description of the malfunction, the dispensed liquids and the complete „Declaration on the Absence of Health Hazards“ to the manufacturer or supplier.

Ask your supplier or the manufacturer for forms. The forms can also be downloaded on www.vitlab.com.

16. Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original parts have been used.

Subject to technical modification without notice.

We will not be held responsible for printing or typographical errors.