



GRIFFIN BEAKER PMP

raised red scale

NEW



- Improved readability
- Enhanced contrast
- Improved durability

GRIFFIN BEAKER

made of PMP & raised red scale

Griffin beakers are used in almost all laboratories, for example to prepare media and for various other applications.

VITLAB's Griffin beakers are a good choice. You get a beaker with **the best durability** combined with **the best readability**.



Better readability due to increased contrast compared to raised scale only.

Better durability compared to the pure printed scale.

Types / Sizes

| Volume ml | Graduation ml | Ø mm | Hight mm | PU | Cat. No. |
|-----------|---------------|------|----------|----|----------|
| 100 | 10 | 55 | 70 | 12 | 60896 |
| 150 | 20 | 66 | 80 | 12 | 60996 |
| 250 | 25 | 77 | 98 | 6 | 61096 |
| 400 | 50 | 87 | 112 | 6 | 61196 |
| 500 | 50 | 94 | 118 | 6 | 61296 |
| 600 | 50 | 100 | 127 | 6 | 61396 |
| 1000 | 100 | 120 | 147 | 6 | 61496 |
| 2000 | 200 | 149 | 187 | 6 | 61596 |
| 3000 | 250 | 170 | 212 | 4 | 61696 |
| 5000 | 500 | 203 | 247 | 4 | 61796 |

VITLAB GmbH
Linus-Pauling-Str. 1
63762 Grossostheim
T +49 6026 977 990
F +49 6026 977 9930

info@vitlab.com
www.vitlab.com



VITLAB is a leading manufacturer of liquid handling instruments and plastic laboratory products worldwide. High class laboratory devices are developed and manufactured in our own production facilities.

Our extensive range of products provides optimal support for your laboratory work in a wide variety of application ranges.

Regardless of whether your work involves volume measurement, titration, dosing, pipetting, sampling or storage: the goal of our development and production is to provide products which facilitate your daily lab tasks and always deliver perfect results.

In addition to the range of products that appears in our catalogue, we also produce plastic products according to customer specification; for example, printed with your company name, your logo or individualised scale.

Our products are premium quality! See for yourself!

BRAND GMBH + CO KG
Wertheim

T +49 9342 8080
info@brand.de
www.brand.de

VACUUBRAND GMBH + CO KG
Wertheim

T +49 9342 8080
info@vacuubrand.com
www.vacuubrand.com

VITLAB GmbH
Grossostheim

T +49 6026 977 990
info@vitlab.com
www.vitlab.com